

Important

Prior to utilizing the Synovasure Online Ordering and Results Portal, ensure that your computer is set-up to allow pop-ups from the portal site

Login to Synovasure Portal (CLICK HERE for portal)

Enter your CD Laboratories Account Number, User ID and Password (case-sensitive) established following first login to the portal. Select 'Log in'. If you are having issues remembering your password, select 'Forgot your password?' to reset password using established security questions.

Log in to your account

[Forgot your password?](#)

Accessing Inbox

On the lefthand toolbar, select 'Inbox' under Laboratory Functions. All results from the last 5 days will be visible in your inbox. Click on patient you intended to view.

User
Inbox (5 days)

Laboratory Functions

Inbox

Patient History

Cumulative Reports

Laboratory Orders

Order Tests

Patient Management

Find Patient

Edit Patient

Add New Patient

Merge Patient

General Information

Activity Logs

Test Reference

Practice Options

User Profile

Print Selected
Display All
Results Per Page: 25
Records 1 - 25 of 39

GO
OK

	Accession Number	Last Name	First Name	DOB	Gender	Collection Date	Physician
-	1234567	MOUSE	MICKEY	01/01/1951	M	07/05/2021	SMITH, M.D., JOHN
	7002 - SYNOVASURE PJI Final						
-	1234579	DOE	JANE	09/05/1928	F	07/02/2021	SMITH, M.D., JOHN
	8260 - CELL COUNT/DIFF, SYNOVIAL FLUID Final						
	7895 - SYNOVASURE MICROBIAL ID PANEL Final						
	7002 - SYNOVASURE PJI Final						
	7120 - SYNOVASURE NEUTROPHIL ELASTASE Final						

Results for requested tests will be presented on the screen. If you would like to download results in a PDF format, click on 'Alternate Report Format' OR select "Print" located above the accession number. A new window will open with the report. If a patient has had prior test results from CD Laboratories, you can view the progression over time by selecting the graph icon next to a result.



Inbox Navigation: Next

Print
Audit
 Reviewed

Accession	Name	DOB	Sex	Req Num	Col Date	DOS
1234567	MOUSE, MICKEY	01/01/1951	M		07/05/2021	07/05/2021
Test Name	Results	Range	Units			
SYNOVASURE PJI	SYNOVASURE ALPHA DEFENSIN PJI	POSITIVE	A			
	ALPHA-DEFENSINS-SF	POSITIVE				
CRP-SF	19.2	H	mg/L			

[Alternate Report Format](#)

If results are older than 5 days...

Select 'Find Patient' under Patient Management on the lefthand toolbar. Fill in the patient's name and select 'Search'.

User

Laboratory Functions

- Inbox
- Patient History
- Cumulative Reports

Laboratory Orders

- Order Tests

Patient Management

- Find Patient**
- Edit Patient
- Add New Patient
- Merge Patient

General Information

- Activity Logs
- Test Reference

Patient Search

Last Name:

First Name:

Patient ID:

Physician:

Accession:

Sex:

Search Date:

OR

Enter a date range (mm/dd/yyyy). -to-

Match Type

	Starts With	Exact
Last Name	<input checked="" type="radio"/>	<input type="radio"/>
First Name	<input checked="" type="radio"/>	<input type="radio"/>
Patient ID	<input type="radio"/>	<input checked="" type="radio"/>
Physician	<input checked="" type="radio"/>	<input type="radio"/>
Accession	<input type="radio"/>	<input checked="" type="radio"/>

Click on the patient ID to view results for requested tests. If patient cannot be found, contact CD Laboratories for assistance.

Search Results

Please select a patient below, or go back and search again.

Practice	Accession	Coll Date	Patient ID	Last Name	First Name	DOB
1234 - ORTHO CLINIC	1234567	07/05/2021		MOUSE	MICKEY	01/01/1951

Need Assistance?

Having trouble navigating the LifePoint Portal for Synovasure results? Can't find your patient? Our team is here to help you and your staff. Give our customer service team a call at **888.981.8378** or email us at **customerservice@cdlaboratories.com**. Be prepared to provide your account number and, if applicable, a specific patient's name and date of birth to help expedite your experience.



888.981.8378



customerservice@cdlaboratories.com